



Where Care Meets Tech



MADE IN BRITAIN

Nurse Call Patient Hand Units

Designed with the user in mind



RNIB endorsed



supporting blind and partially sighted people

Nurse Call Patient Hand Units

Designed with the user in mind

Iconic design

A well designed hand unit is essential for patient well being. Static Systems' patient hand units are amongst an elite group of products that can claim to be iconic. The design set new standards and has become the benchmark product in its field.

Patient-centric design

The design of our patient hand units was achieved following consultation with user groups and professional bodies, including the Royal National Institute of the Blind (RNIB). Units have an easy-grip profile, reassuringly tactile buttons and distinguishing colours. All adding up to a thoroughly researched design inspired by a simple aim – to bring ease of use and peace of mind to patients.

Round-the-clock protection

To assist in infection control measures, enhanced protection against bacteria is provided by a silver-based antimicrobial agent incorporated during the manufacturing process. Units are sealed to IP67 and can be cleaned with a sterilizing wipe to eliminate bio-hazards and reduce cross-infection risks between patients.

Extended life cycle

Units are designed to be robust and essentially maintenance-free, comprising a minimal number of component parts. Damaged cables can be easily replaced.

Outstanding features...

Stylishly contoured to fit comfortably in the hand, the carefully considered design provides a number of features that enable the patient to quickly identify their needs and make confident choices, whilst the robust construction ensures suitability for use within the hospital environment.

Acoustically designed for high quality speech


Sealed microphone and speaker

Distinguishing colours for patient ease of use, as recommended by the RNIB

Tactile buttons for user reassurance

Push buttons control the patient's lighting environment

Radio and TV functions provided for in-bed entertainment



Patented quick-release plug prevents injury to the patient and damage to equipment

Cable strain relief to minimise potential cable damage

Replaceable cable extends product life

Linen clip prevents the hand unit from slipping out of reach



The complete family of advanced nurse call solutions



Ultima™ / Full end-to-end, open protocol IP solution



Aspire smart_{sync}® / Wireless solution



Codemlon™ / Addressable loop wired solution



Where Care Meets Tech

www.staticsystems.co.uk

Static Systems Group Limited
Heath Mill Road, Wombourne, Staffordshire, WV5 8AN, UK.
T: +44 (0) 1902 895 551

In line with our company policy of continuous product development, we reserve the right to change design and improve specification without prior notice.

Focus:008 (02/2021)



A **Halma** company